

PP Non-woven Face Mask

3-ply, earloop, blue, Type II R

Product Details

- 3 plies
- latex-free elastic ear loops
- integrated, adjustable nose bar
- suitable for person wearing glasses
- blue

Material

- inner layer: PP spunbond non-woven, white
- filter medium: PP meltblown non-woven, white
- outer layer: PP spunbond non-woven, blue
- skin-friendly, non-irritating, light

Material Weight

- inner layer $\approx$ 25 gsm ($\pm$ 2 gsm)
- filter medium $\approx$ 28 gsm ($\pm$ 2 gsm)
- outer layer $\approx$ 30 gsm ($\pm$ 2 gsm)

Filtration Performance (BFE)

- $\geq$ 99%

Norms & Guidelines

- CE according to Medical Device Directive 93/42/EEC + 2007/47/EC – Class I product
- Disposable product
- Latex-free
- Non sterile
- EN 14683:2019 + AC:2019
- ISO 13485:2016

Packaging

50 pcs/dispenser, 2000 pcs/carton

