

Specification Sheet


PE coated isolation Gown - Blue

Product Details

- size 115x137cm
- Blue
- PE complete coating
- knitted cuffs
- ties in neck and waist

Material

- spunbond
polypropylene+polyethylene
- solid, light
- moisture-resistant
- ≈ 30 gsm (±3 gsm)



Norms & Guidelines

- CE according to Medical Device Directive 93/42/EEC + 2007/47/EC – Class I product
- Disposable product
- Latex-free
- Non-sterile
- ISO 13485:2016

Packaging

10 pcs/polybag, 100 pcs/carton

Code	Color	Export per carton	Gross Weight per carton	Net Weight per Carton	Carton Dimension
30052	blue	100 pcs	7,2 kg	6,5 kg	52 x 28 x 33 cm



Subject to modification without prior notice

Jul., 2020