

Specification Sheet


PP Clip cap - single elastic band

Product Details

- 53cm
- white, blue, green, red, yellow, coral
- single elastic band

Material

- spunbond polypropylene
- breathable, light
- ≈ 10 gsm (± 1 gsm)



Norms & Guidelines

- CE according to Medical Device Directive 93/42/EEC + 2007/47/EC – Class I product
- Disposable product
- Latex-free
- Non-sterile
- ISO 13485:2016

Packaging

100 pcs/polybag, 1000 pcs/carton

Code	Color	Export per carton	Gross Weight per carton	Net Weight per Carton	Carton Dimension
40131	white	1.000 pcs	3,1 kg	2,7 kg	34 x 25 x 24 cm
40132	blue	1.000 pcs	3,1 kg	2,7 kg	34 x 25 x 24 cm
40133	green	1.000 pcs	3,1 kg	2,7 kg	34 x 25 x 24 cm
40134	red	1.000 pcs	3,1 kg	2,7 kg	34 x 25 x 24 cm
40135	yellow	1.000 pcs	3,1 kg	2,7 kg	34 x 25 x 24 cm
40136	coral	1.000 pcs	3,1 kg	2,7 kg	34 x 25 x 24 cm



Subject to modification without prior notice

Jan., 2021

Dispo-Dent B.V. E-mail: info@dispo-dent.nl Tel: 0321-318140 Mobile: 0651418777