


Specification Sheet

HDPE Sleeve cover - Blue

Product Details

- size 20(+2cm)x40 cm
- ≈0.016mm, 2.8g/pc
- blue
- smooth surface
- machine-made



40 x 20cm, 0.016mm

Material

- high density polyethylene
- moisture-resistant
- skin-friendly, non-irritating, light

Norms & Guidelines

- CE according to Medical Device Directive 93/42/EEC + 2007/47/EC – Class I product
- Disposable product
- Latex-free
- Non sterile
- ISO 13485:2016
- Suitable for contact with food

Packaging

100 pcs/polybag, 2000 pcs/carton

Code	Color	Export per carton	Gross Weight per carton	Net Weight per Carton	Carton Dimension
60116	blue	2000 pcs	6,0 kg	5,6 kg	40 x 22 x 32 cm



Subject to modification without prior notice

May, 2020