

Specification Sheet


HDPE Glove W/paper card- Embossed, Transparent

Product Details

- 250mm x 280mm (+ 20mm)
- ≈0,013mm, 0,85g/pc
- transparent
- embossed surface
- ambidextrous
- with paper card

Material

- high density polyethylene
- moisture-resistant
- skin-friendly, non-irritating, light



Norms & Guidelines

- CE according to Medical Device Directive 93/42/EEC + 2007/47/EC – Class I product
- Disposable product
- Latex-free
- Non sterile
- ISO 13485:2016
- Suitable for contact with food

Packaging

50 pcs/paper card, 500 pcs/polybag, 5,000 pcs/carton

Code	Color	Export per carton	Gross Weight per carton	Net Weight per Carton	Carton Dimension
60224	transparent	5.000 pcs	5,7 kg	4,3 kg	36 x 22,5 x 20 cm



Subject to modification without prior notice

Jul., 2020