


Specification Sheet

LDPE Apron on block - Blue

Product Details

- size 90x160 (+3) cm
- ≈0,05mm, 50g/pc
- on block
- blue
- embossed surface
- neck loop
- ties in waist



Material

- low density polyethylene
- moisture-resistant
- skin-friendly, non-irritating, light

Norms & Guidelines

- CE according to Medical Device Directive 93/42/EEC + 2007/47/EC – Class I product
- Disposable product
- Latex-free
- Non sterile
- ISO 13485:2016
- Suitable for contact with food

Packaging

25 pcs/polybag, 250 pcs/carton

Code	Color	Export per carton	Gross Weight per carton	Net Weight per Carton	Carton Dimension
60342	blue	500 pcs	13 kg	12,5 kg	47,5 x 30 x 13,5 cm



Subject to modification without prior notice

Jul., 2020